

# Fulflo® Poly-Mate™ Plus Filter Cartridges

High surface area and high efficiency  
all-polypropylene pleated cartridges

Fulflo® Poly-Mate™ Plus Cartridges, made of pleated polypropylene microfiber, provide high efficiency and high purity filtration. The high efficiency of the Poly-Mate Plus line makes it an ideal membrane pre-filter or cost-effective alternative to membrane cartridges in a wide range of applications.

Poly-Mate Plus Pleated Cartridges are available in the following pore sizes (nominal rating at 90%): 0.25µm, 0.45µm, 0.8µm, 2.0µm, 3.0µm, 5.0µm, 30.0µm, 50.0µm, 100.0µm.



## Contact Information

Parker Hannifin Corporation  
**Bioscience Division - N.A.**  
2340 Eastman Avenue  
Oxnard, CA 93030

phone +1 805 604 3583  
bioscience.na@parker.com

[www.parker.com/bioscience](http://www.parker.com/bioscience)

## Benefits

- All-polypropylene media and construction meet a broad range of performance requirements
- One-piece integral construction is 100% bonded for maximum cartridge integrity
- High surface area design provides superior flow rates and extended service life
- All materials of construction are FDA listed as acceptable for potable and edible liquid contact according to CFR Title 21
- Fixed pore construction provides ultimate particle retention
- Major end seal options are available to fit most standard vessels
- Poly-Mate™ Plus cartridges are non-fiber releasing and ensure consistent quality filtration performance
- ISO 9001 registered company

## Applications

- DI Water
- Process Water
- Magnetic Media
- Plating Chemicals
- Membrane Pre-filter



ENGINEERING YOUR SUCCESS.

# Fulflo® Poly-Mate™ Plus Filter Cartridges

## SPECIFICATIONS

### Materials of Construction

#### Filter Media

- Melt blown polypropylene microfiber

#### Media Support Layers

- Non-woven or mesh polypropylene

#### Core

- Heavy wall high strength polypropylene

#### Media Support Cage and Thermally

#### Welded End Caps

- Molded polypropylene

#### Seal Materials

- Buna-N, EPR, Silicone, Viton®, PFA Encapsulated Viton®

### Dimensions:

#### Cartridge Outside Diameter:

- 2 1/16 in.

#### Cartridge Inside Diameter:

- DOE: 1 1/16 in.
- SOE: 1 3/32 in.

### Maximum Recommended Operating Conditions:

Temperature - 200°F (93°C)

Temperature @ 35psid - 160°F (71°C)

Change Out ΔP - 35psi (2.4bar)

ΔP @ Ambient 70°F (21°C) - 70psi (4.8bar)

ΔP @ 200°F (93°C) - 20psi (1.4bar)

Flow Rate -10gpm (38 lpm) per 10 in. length

### Product Safety:

- All components FDA listed per CFR, Title 21
- Non-fiber releasing per FDA Part 210.3B (5) and (6)
- Non-photosensitive

### Filtration Ratings:

90% at 0.25, 0.45, 0.8, 2, 3, 5, 10, 30, 50 and 100 micrometer pore sizes

### Performance Attributes

#### Flow Rate and Pressure Drop Formulas

$$\text{Flow Rate (gpm)} = \frac{\text{Clean } \Delta P \times \text{Length Factor}}{\text{Viscosity} \times \text{Flow Factor}}$$

$$\text{Clean } \Delta P = \frac{\text{Flow Rate} \times \text{Viscosity} \times \text{Flow Factor}}{\text{Length Factor}}$$

#### Notes:

- Clean ΔP is psi differential at start.
- Viscosity is centistokes. Use Conversion Tables for other units.
- Flow Factor is psid/gpm at 1 cks for 10 in. (or single).
- Length Factors convert flow or ΔP from 10 in. (single length) to required cartridge length.

### Liquid Particle Retention Ratings (μm) @ Removal Efficiency of:

Cart.	β=1000 99.9%	β=100 99%	β=50 98%	β=20 95%	β=10 90%
PMP002	2.2	1.6	0.90	0.45	0.30
PMP004	3.1	2.9	1.4	0.75	0.45
PMP008	9.2	8.0	3.2	1.5	0.8
PMP020	11.0	9.5	8.6	3.1	1.7
PMP030	12.0	11.0	6.1	4.6	3.0
PMP050	14.0	12.0	10.6	8.4	5.0
PMP100	21.0	17.0	15.0	12.0	10.0
PMP300	52.0	44.0	35.0	24.0	15.0
PMP500	71.0	68.0	62.0	56.0	50.0
PMP1000	138.0	126.0	117.0	109.0	100.0

### Poly-Mate Plus Flow Factors (psid/gpm @ 1 cks)

Rating (μm)	Flow Factor
0.25	0.0900
0.45	0.0530
0.8	0.0290
2	0.0068
3	0.0060
5	0.0048
10	0.0040
30	0.0030
50	0.0025
100	0.0020

### Poly-Mate Plus Length Factors

In.	Factor
4	0.4
10	1.0
20	2.0
30	3.0
40	4.0

## Ordering Information

Cartridge Code		Pore Size		Nominal Length		Support Construction		Seal Material		End Cap Configuration		Special Options	
CODE	MATERIAL	CODE	MICRON	CODE	INCHES	CODE	MATERIAL	CODE	MATERIAL	CODE	DESCRIPTION	CODE	DESCRIPTION
PMP	Poly-Mate Plus	002	0.25	4	4	A	Natural Polypropylene (All support components)	E	EPR	AR	020 O-ring/Recessed cap	No Symbol	No Option
		004	0.45	10	9 1/16			N	Buna-N	DO	Double open end (DOE)	B	Bubble-point test
		008	0.8	20	19 1/16			S	Silicone	DX	Double open end/extended core	R	DI water rinse (5 minutes)
		020	2	30	30 1/8			T	PFA/Viton® (SOE)	LL	120/120 (Filterite LMO & Nuclepore Polymeric Vessels)**	Z6	Individual Poly bag only
		030	3	40	40			V	Viton®	LR	120 O-ring/Recessed (Nuclepore)**		
		050	5							PR	213 O-ring/Recessed cap (Ametek® & Parker LT Polymeric Vessels) **		
		100	10							TC	222 O-ring/Flat		
		300	30							TF	222 O-ring/Fin		
		500	50							SC	226 O-ring/Flat		
		1000	100							SF	226 O-ring/Fin		

\*\*Available only in 9 3/16" (-9) and 19 3/16" (-10) lengths

Specifications are subject to change without notification.  
For User Responsibility Statement, see www.parker.com/safety

© 2017 Parker-Hannifin Corporation  
Bioscience Division - North America  
All Rights Reserved



Fulflo is a registered trademark of Parker-Hannifin Corporation  
Viton is a registered trademark of E.I. DuPont de Nemours & Co., Inc.  
Ametek is a registered trademark of Ametek, Inc.

DS\_IP\_Poly-Mate Plus Rev. C