

Scaleable SINGLE-USE filtration solutions



Effective Filtration Area (EFA)	Syringe Filters		E Size		B Size		A Size		K Size		10" Size	
	cm ²	in ²	m ²	ft ²	m ²	ft ²	m ²	ft ²	m ²	ft ²	m ²	ft ²
TETPOR AIR	14.50	2.25	0.06	0.64	0.12	1.29	0.25	2.69	0.36	3.87	0.77	8.28
HIGH FLOW TETPOR II	14.50	2.25	-	-	-	-	-	-	-	-	-	-
PROCLEAR GF	14.50	2.25	0.05	0.6	0.10	1.07	0.20	2.2	0.27	2.9	0.56	6.0
PROCLEAR GP	14.50	2.25	0.03	0.3	0.06	0.6	0.12	1.3	0.16	1.7	0.34	3.7
PROCLEAR PP upto	14.50	2.25	0.07	0.75	0.14	1.50	0.28	3.01	0.37	3.98	0.79	8.5
PROPOR BR	14.50	2.25	0.05	0.53	0.10	1.07	0.20	2.15	0.26	2.79	0.55	5.92
PROPOR SG	14.50	2.25	0.05	0.53	0.10	1.07	0.20	2.15	0.26	2.79	0.55	5.92
PROPOR HC	14.50	2.25	0.05	0.53	0.10	1.07	0.20	2.15	0.26	2.79	0.55	5.92
PROPOR LR	14.50	2.25	0.05	0.53	0.10	1.07	0.20	2.15	0.26	2.79	0.55	5.92
PROPOR MR	14.50	2.25	0.05	0.53	0.10	1.03	0.19	2.09	0.24	2.58	0.50	5.38

